


**600026 HOME CARE BED WOODEN SIDE
USER MANUAL**



600026



MOTORISED HOME CARE PATIENT BED

Contents

1. Introduction	3
2. PRODUCT DESCRIPTION	3
3. SAFETY AND WARNING	3
 3.1. Safety and Warning Instructions	3
4. SYMBOLS AND LABELS ON THE PRODUCT	4
5. AREA AND CONDITIONS OF USE	6
6. APPLIED STANDARDS AND REGULATIONS	6
 6.1. Compliance with Standards	6
 6.2. Gap Distances	7
 6.3. Compression Distances	7
 6.4. Foot Impingement Distance	7
 6.5. Balance and Safe Working Load	8
 6.6. Distance between Mattress and Side Railing	9
 6.7. Fall Protection	9
 6.8. Side Railing Strength	9
 6.9. Angular Movements	10
7. SYMBOLS RELATED TO CONTROL TOOLS AND PERFORMANCE	10
8. THE PURPOSE AND INSTALLATION OF THE PATIENT BED	11
 8.1. Transfer	11
 8.2. Installation and Assembly	12
9. TRANSPORT OF THE PATIENT	14
 9.1 Considerations for Patient Transport	14
10. OPERATION	14
 10.1. Initial Start-up	14
11. ELECTRICAL SYSTEM AND CONTROL ELEMENTS	14
11.1 Electrical System	14
 11.2. Movement Functions of the Patient Bed	15
 11.2.1. Hand Control	16
 11.2.2. Side Railing	16
 11.2.3. Brake System	16

11.2.4 Cross Brake System	16
11.2.5. Foot Platform Angular Movement (Ratchet Movement)	16
12. TECHNICAL DIMENSIONS	17
MAINTENANCE, REPAIR AND CLEANING INSTRUCTIONS.....	18
13. Maintenance and Repair Periods	18
13.1 Monthly Maintenance	18
13.2 Annual Maintenance	18
14. Spare Parts	18
15. Cleaning / Disinfecting Instruction.....	19
15.1 Cleaning	19
15.2. Disinfection.....	19
16. TROUBLESHOOTING	19
17. SAFE STORAGE OF THE PATIENT BED.....	20
18. ENVIRONMENT.....	20
18.1 Environmental Protection	20
19. WARRANTY.....	20
20. CONTACT	20

1. INTRODUCTION

The use of the electric patient bed used for adults, cleaning, technical information, production design in accordance with standards and regulations, and the use of the product is not affected by any risk and product performance.

2. PRODUCT DESCRIPTION

NHB 7000 Models is an electric patient bed with 4th motorised back, foot, height and Trendelenburg movement functions designed for safe treatment of the patient. The bed can take Fowler/Vascular positions.

600026 Models are designed based on maximum safety and mobility. 600026 models, which are designed by taking into account the physical disabilities and comfort of the patient, maximise the comfort of patients in the high risk group in every sense.

600026 Model electric patient beds are designed in accordance with international standards. Thus, high performance is provided. The patient's fall, jamming, etc. The design has been carried out by minimising the negativeities that may occur against risks.

3. SAFETY AND WARNING


3.1. Safety and Warning Instructions



- Follow the instructions carefully and, if necessary, keep a copy of the instructions close to the bed.
- To avoid customer damage and to obtain the best performance from the product, this user manual must be read before use.
- Use the correct mains supply for the bed
- Make sure that the bed is operated by qualified personnel
- To eliminate the risk of electric shock, existing electrical medical equipment may only be connected to a supply network with protective earthing.
- If the power supply cables are not handled properly, there is a danger of tangling, cutting or other mechanical damage.
- While the patient bed functions are in motion, the cables may be pinched, so do not perform any functionality on the bed by getting technical support in case of sagging cables.
- Do not interfere with the bed while the functions of the patient bed are in motion, it poses a danger due to limb entrapment of fingers, head, etc. through existing gaps.
- Do not use any other mattress other than the 860x1950x12 mm 32 density mattress used by the manufacturer. Dangers such as falling and jamming occur due to mattress dimensions.
- The distance between the top surface of the side rail and the top surface of the mattress should be at least **220mm** when the side rail is locked.
- To prevent damage to the POWER SUPPLY CORD, wrap the POWER SUPPLY CORD around the cable winder mechanism that keeps the POWER SUPPLY CORD away from any moving parts or mechanisms of the PATIENT BED when the PATIENT BED is being used, transported or not in use.
- Lift the patient bed max. to a height of **750 mm** from the floor.
- Do not apply a load of more than **175 kg** to the patient bed.
- Only original spare parts and damaged parts should be replaced immediately and parts should be replaced by qualified personnel
- Do not overload the bed beyond its carrying capacity.
- If you notice that the bed does not fulfil its functions, do not use the bed and inform the technical service.
- Inform the patient and other users about the operating functions.
- The patient bed should be used indoors (wards, observation rooms, etc.).
- In case of any electrical and mechanical problems in the patient bed, intervention to the bed should be carried out by trained personnel. If necessary, the manufacturer's technical service unit should be informed and support should be obtained.
- During the delivery of the patient bed, training should be requested from the manufacturer for product use and intervention in case of any malfunction.
- Detailed technical information about the patient bed should be requested if necessary.
- The safe operating weight must never be exceeded. If the safe operating weight must be exceeded, the surface on which the patient lies must be kept in its lowest position and the functions of the bed must not be used.
- The 2nd person should not sit on the bed while the patient is lying down.
- To reduce the risk of injury from falling off the bed when the patient is left alone, the bed should be raised when it is at the bottom.
- No other non-production independent mechanism should be connected to the patient bed.

- In case of any malfunction, it should not intervene except for the authorised technical service.
 - The bed must never be used by more than one patient at the same time.
 - When setting the bed in the desired position, the condition of the patient and the bed environment should be taken into account.
 - Always unplug the power cord before cleaning or replacing any part of the bed.
 - Deliver the defective motor, plastic, etc. materials used in the patient bed to the manufacturer or licensed waste company in terms of environmental protection.
 - Deliver the packaging material such as cardboard, nylon, etc. used in product packaging to the licensed organisation or manufacturer company for environmental protection.
 - The patient or carer is definitely in danger if
 - ❖ When the power cable is damaged
 - ❖ When moving the bed from one place to another when the floor and bed surroundings are unreliable.
 - ❖ Improper maintenance (e.g. automatic washing or washing with pressurised water)
 - ❖ When the safe working weight is exceeded.
 - Only spare parts supplied by an authorised workshop should be used. If spare parts supplied by other suppliers are used, the bed frame manufacturer does not accept any responsibility for any damage, loss or injury.
 - The patient bed should be wiped with a damp cloth with a cleaning and disinfection solution with water at room temperature.
-


4. SYMBOLS AND LABELS ON THE PRODUCT

	Protection against accidents due to electricity type B		Bed load 220 kg
	CE marking		IEC 60417-5019 / Earthing
	Accompanying documents		Suitable for indoor use.
	Read the user manual		General Warning Sign
	Equipotentiality		CPR
	Safe payload: 220 kg		Entrapment hazard
	Maximum patient load: 175 kg		Label used where removable parts exceed 20 kg
	Bed dimensional hazard, read the operating instructions, siderail functional hazard and caution symbol S		Company contact information, barcode number, product model, serial number, production date, label with medical device class



If you see that any of the labels on the product is missing, please request it from the manufacturer.

WARNING



1	FootEndPanel
2	Foot platform
3	Side Railing
4	Back Platform
5	HeadEnd Panel
6	Side railing opening and closing latch
7	Ø125 Wheel
8	Brake Pedal

5. AREA AND CONDITIONS OF USE

Patientbeds; It is designed for the treatment of adults according to usage areas such as hospitals, homes and closed environments where medical interventions are performed. Using patientbeds outside of the instructions for use can cause serious injuries and poses a danger.

The patientbed includes all parts of the patientbed accessible to the patient, even if the part to which the patientbed is applied is under the mattress support platform.



Use the patientbed in the areas approved by the manufacturer and according to the production purpose.

WARNING

The use of the patientbed may pose a hazard depending on temperature and humidity. Use the bed in the following environments. Temperature Range -10°C to +40°C

- Relative humidity between 30 and 75
- The bed is intended for use inside rooms for medical purposes. Electrical installations must therefore meet compliance with local connections.

6. APPLIED STANDARDS AND REGULATIONS


The patientbed complies with the following standards and directives.

- TS EN ISO 14971
- 93/42/EEC Replaced by MDR 2017/745 EU
- AT DECLARATION OF CONFORMITY


6.1. Compliance with Standards

Our patientbed, which we produce within the framework of international standards, keeps high performance, ergonomic and safety at the highest level.

Examples of patient entrapment in the patientbed and the measures taken in this direction;




6.2. GapDistances



REGION	VALUE
A	$\leq 120 \text{ mm}$
B	$\leq 60 \text{ mm}$
C	$\leq 60\text{mm Or} > 318 \text{ mm}$

6.3. JammingDistances




6.4. FootImpingementDistance

Description

For the zone where dimension "b" is less than or equal to 130 mm, dimension "a" is always greater than or equal to 120 mm.

. Description

For the zone where dimension "b" is between 130 mm and 180 mm, dimension "c" is always greater than or equal to 50 mm




6.5. Balance and Safe Working Load

Static loads of **175 kg** applied to the patient bed are applied to the areas required by the standard as shown in the figure below and the stability of the patient bed is tested.


Safe carrying capacity of the patient bed is **175 kg**.

Static carrying capacity is **350 kg**.



The SAFE OPERATING LOAD of the BED LIFT must be at least **1750 N**. This load is accepted as the sum of the following minimum loads.

- ✳ 1350 N, corresponding to a mass of approximately **135 kg** for the PATIENT,
- ✳ 200 N, corresponding to a mass of approximately **20 kg** for the mattress,
- ✳ 150 N, corresponding to a mass of approximately **15 kg** for the SAFE WORKING LOAD supported by these ACCESSORIES, excluding ACCESSORIES and PATIENT mass.
- ✳ 500 N, corresponding to a mass of approximately **50 kg** for those sections of the PATIENT BED designed to be lifted by the BED LIFT,



REGION	VALUE
A	45% OF THE SAFE WORKING LOAD OF THE BACK SECTION
B	25 % OF THE SEATING SECTION SAFE WORKING LOAD
C	30% OF THE LEG SECTION SAFE WORKING LOAD
SWL	SAFE WORKING LOAD


6.6. Distance Between Mattress and Side Railing

Positive results are obtained in cases where the distance between the mattress support platform and the side railing does not exceed 50% of the **120 mm - 60 mm** conical gauge.

A force of **250 N** is applied to the **60 mm** cylindrical end of the conical tool in the most unfavourable direction. The large end of the conical tool must not sink below the mattress surface by more than **50 %** or more than **120 mm** in diameter.

6.7. Fall Protection


Patient bed side rails and mattress top surface are designed with minimum height rules. In this way, the risk of the patient falling is minimised by making a risk assessment and the design is carried out in accordance with the standards.




REGION	VALUE
1	BEDSIDE PANEL
2	FOOTPAD PANEL
3	SHEET
G	$\geq 220 \text{ mm}$

6.8. Side Railing Strength


Forces are applied to the side rails in accordance with EN 60601-2-52 standard.






REGION	APPLIED FORCE
A	500 N
B	500 N
C	750 N
D	750 N
E	500 N
F	500 N



6.9. Angular Movements



REGION	ANGLE
A (Angle between the backrest and the level plane)	0° - 75°
B (The joint point of the line drawn between the backrest section and the section/sitting section where the lower part of the leg is placed and the section/sitting section where the upper part of the leg is)	Min 90°
C The angle B between the horizontal plane and the line drawn between the point of articulation of the backrest section/seating section and the point of articulation of the upper leg section/lower leg section is at least between 0° and 12°.	0° - 32°
D (Trendelenburg angle) must be adjustable	Min 12°

7. SYMBOLS RELATED TO CONTROL TOOLS AND PERFORMANCE



	Upwardfunction of thelegsection 		Downwardfunction of thelegsection 
	Downwardstraightpositionmovementfun ction 		Upwardstraightpositionmovementfuncti on 
	Reversetrendelenburg 		Straighttrendelenburg 

8. THE PURPOSE AND INSTALLATION OF THE PATIENT BED

8.1. Transfer

Considerations for a safe transport;

- Make sure that no cables are disconnected during the transfer.
- Make sure that the power cable is wrapped around the hook at the bedside of the patient bed
- Make sure that the wheels are locked when loading and unloading
- Transport the bed on a suitable surface.
- Check that the brakes are on during transport.
- Remove any accessories that may fall off during movement.



Do not forget to disconnect the power cord from the socket before transporting the patient bed.

ATTENTION

8.2. Installation and Assembly

Adjust the bed as follows

- The following rules must be followed during the installation of the bed.
- Please refer to the scope of delivery and bed variants.
- Attach the accessories accordingly.
- Make sure that the ground is suitable during the installation process.
- Make sure that the connection cables are inserted into the correct sockets.
- Check the earthing cables.
- Deliver the dismantled packaging to the necessary places for environmental health.


Product Folding:



Please act in accordance with the instructions in order not to cause limb entrapment during product folding.

ATTENTION

1. Step:



Remove the handle by turning it anti-clockwise.

2. Step:




After removing the handle, remove the siderails.

3. Step:



After performing step 2, remove the handles at the marked points.

4. Step:




During the folding process, hold on to the area marked in red, there may be a risk of limb entrapment due to moving parts.

Fold from the empty brake part by activating the brake system of any of the wheels from the parts where the head or foot panels are located.

5. Step:



After folding, fix any of the removed handles to the marked place by turning clockwise.



Reassemble the bed by following the steps above in reverse order.



NOTE: If any problem is encountered during the installation of the patient bed, the intervention should be carried out by PEDIA PALS technical service personnel.

Caution Risk of injury

During installation, install moving accessories last, as there is a risk of injury from falling.

When installing the electrical system, check the cables and install them. Any deformation of the cables poses a risk of electric shock.

Ensure that the main supply is disconnected during bed assembly. Ensure that the castors are locked before assembly.

Installation of the patient bed should be carried out by technical service personnel or trained hospital personnel.

9. TRANSPORT OF THE PATIENT

9.1 Considerations for Patient Transport

Make sure that the patient bed is in the appropriate position.

- Make sure that no cables are connected during the transfer.
 - Make sure that the power cable is wrapped around the hook at the bedside of the patient bed
 - Make sure that the wheels are locked when loading and unloading
 - Transport the bed on a suitable surface.
 - Check that the brakes are on during transport.
 - Remove any accessories that may fall off during movement.
 - Make sure the patient is in the centre of the bed.
 - Make sure that the side rails are above.
 - Keep the bed at the lowest level.
-
-

10. OPERATION

10.1. Initial Start-up

Prepare the bed as follows;

- Check the electrical connection points (socket inputs) before inserting the plug of the mattress into the socket outlet.
 - Check the connections of the grounding cables.
 - Check the insulation of the connection cables. Check for any deformation or crushing.
 - Check the wheels and also check the function of the brake system wheels.
 - Plug the power cable into the socket for power supply. Make sure the power is switched on
 - Operate and check the back, foot, height and other function motors.
 - Check the control
 - Check the function of the side rails.
-

11. ELECTRICAL SYSTEM AND CONTROL ELEMENTS

11.1 Electrical System

The patient bed has a control box, foot motor, back motor, calf movement, height motor and hand control. As a working principle, the system consists of the control box receiving the data and performing the functions of the functions with the command process given from the hand control. Motors Below is the connection diagram of the system. Safety rules must be followed during connection. Read the instructions carefully in this regard. Read the electrical information on the product label. The components and structure of 600026 models are designed in accordance with safety rules. The models have IPX4 / IPX6 degree of protection.



WARNING: Electrical interventions should be carried out by authorised and trained personnel and support should be obtained from the manufacturer when necessary.


CAUTION: Do not touch the control box, motor and cables when the bed is in operation.

CAUTION: If the patient bed is connected to an off-grid power supply, check the periodic checks of the power supply and see instructions.

CAUTION: Electrical intervention must be carried out by trained personnel. Accidents due to electric shock may occur.

Electrical System Connection Diagram

ÖGE	PARÇA NO	ADET
1	CA4 KONTROL KUTUSU J11332	1
2	009914948-1700-B	4
3	J20801 LINAK SIRT MOTORU	1
4	J01411 LINAK YÜKSEKLİK	2
5	J20800 LINAK AYAK	1
6	J14051 HB850200130004	1
7	SML912263	1



ELECTRICAL PROPERTIES

COMPONENT	MANUFACTURER / BRAND	TECHNICAL INFO
FOOT MOTOR	LINAK	IMPORTAL / 6000 N / DC / IPX4 / IPX6 / 24V / Max. 5.0 A
BACK MOTOR	LINAK	IMPORTAL / 3500 N 4500 N / DC / IPX4 / IPX6 / 24 V / Max. 3.5 A
HEIGHT MOTOR	LINAK	IMPORTAL / 6000 N / DC / IPX4 / IPX6 24V / Max. 5.0 A
CONTROL BOX	LINAK	IMPORT / 100-240 V ~ / 50/60 Hz / Max. 5 A / IPX6
HAND CONTROL	LINAK	IPX4 IPX6

11.2. Movement Functions of the Patient Bed

The patient bed performs back angular movement, foot angular movement, height and trendelenburg movements with the help of an electric motor by giving a command from the hand control. In the foot section, the foot section is manually moved from the angular position to the parallel position by means of a 5-stage ratchet. The CPR arm on the patient bed also fulfills its function manually.

The siderails on the patient bed fulfill their functions manually with the help of the lock mechanism.



Dangers of injury when performing the movement function of the patient bed

ATTENTION

Make sure that there are no body parts between the patient bed platform. When adjusting the positions of the patient bed, stay away from pinch points. Do not move the patient bed beyond the safe load.


Model-dependent control elements;

- Wheel

11.2.1. Hand Control



The bed is designed for movement control. It is connected with a flexible cable and can be easily used by the carer and the patient in any position. Detailed function movements are given in Article 7.

11.2.2. Side Railing



Pull the siderail latches inwards and press down slightly on the siderail.



WHEN THE SIDE RAILING IS LIFTED DIRECTLY UPWARDS, IT COMES TO THE CLOSED FUNCTION.



11.2.3. Brake System


Brake system: Designed using the pedals on the wheel.

11.2.4 Brake Pedal





WHEN THE BRAKING PEDAL ON THESE WHEELS IS PRESSED DOWNTOWARDS, THE BRAKE IS ACTIVATED. WHEN THE PEDAL IS MOVED UPWARDS, THE BRAKES BECOME INACTIVE.

11.2.5. Foot Platform Angular Movement (Ratchet Movement)




FOOT PLATFORM ANGULAR MOVEMENT IS REALISED GRADUALLY WITH THE HELP OF CIRCI. WHEN YOU LIFT THE PLATFORM UPWARDS, ANGULAR MOVEMENT OCCURS GRADUALLY BY MOVING THE CIRCI TO THE TOP POSITION. WHEN YOU MOVE THE PLATFORM DOWNTOWARDS, THE CIRCI IS MOVED TO THE TOP POSITION AND THEN THE PLATFORM IS RELEASED FREELY.

12. TECHNICAL DIMENSIONS



WHEEL DIAMETER IS 125 MM



600026 MOTORIESD HOME CARE PATIENT BED TECHNICAL SPECIFICATIONS

		Measure	Unit	
A	BedFrameLength (OutsidetoOutside)	2241	mm	
B	TiltSurfaceHeightAboveGround	Max. 750 - Min. 350	mm	
C	BedFrameWidth (OutsidetoOutside)	1022	mm	
D	Height of IV Pole	min.	1100	mm
		max.	1250	mm
α	BackAngle	0-75	$^{\circ}$	
β	CalfAngle	0-32	$^{\circ}$	
θ	FootAngle	0-12	$^{\circ}$	
	SafeLoadCapacity	175	Kg	
	StaticCarryingCapacity	350	Kg	
	Wheel Diameter	Ø125	mm	
	TrendelenburgAngle	12	$^{\circ}$	
	Mattress	32	d	

MAINTENANCE, REPAIR AND CLEANING INSTRUCTIONS

13. Maintenance and Repair Periods



ATTENTION

Injury may result from improper maintenance.

- ✓ If the maintenance supervisor is unsure, he/she should seek support from the manufacturer.
- ✓ Take necessary safety precautions before maintenance.
- ✓ Do not use spare parts not recommended by the manufacturer.
- ✓ Do not perform maintenance while the bed is in functional movement.

NOTE: maintenance records of patient beds should not be kept and these intervals should be made at certain times.

13.1 Monthly Maintenance

- Check the moving parts of the patient bed (back, height, trendelenburg movement, etc.).
- Check the joining elements (bolts, nuts, etc.)
- Check the accessory slots for any wear and tear.
- Check the wheels and their function.
- Check the brake system.
- Check the side rail movement functions.

13.2 Annual Maintenance

- Check the joints of the movement functions.
- Check the bed frame joints.
- Check the brake system bolts and joints.
- Check the functions of the hydraulic parts and check for any leakage or leakage.
- Check the wheel function and check the wheel mounting bolts.
- Check the side railing connection points and fasteners.



WARNING

- ✓ Malfunctions caused by errors in use and spare parts demand are excluded from the scope of warranty.
- ✓ Maintenance should be carried out by trained hospital staff, if any problems are encountered, our company's technical service unit should be informed.
- ✓ Check all bolts and tighten if necessary
- ✓ Replace worn accessories
- ✓ Do not replace worn materials with faulty materials.

CAUTION: The patient bed may be damaged by improper maintenance.

If the authorised personnel of the hospital are unsure, seek advice from PEDIA PALS. Maintenance should only be carried out by authorised, trained personnel.

Receiving spare parts requests and information

14. Spare Parts

Malfunctions caused by errors in use and spare parts demand are excluded from the scope of warranty.

Maintenance should be carried out by trained hospital staff, if any problems are encountered, our company's technical service unit should be informed.

Receiving spare parts requests and information

- ✓ PEDIA PALS Technical service
- ✓ service@pediapals.com
- ✓ 1-888-733-4272

15. Cleaning / Disinfecting Instruction

15.1 Cleaning

- Uses suitable detergents for cleaning. The VOC values of these detergents should be at appropriate values that will not harm the environment and people.
- Do not use abrasive powders, steel wool, steel wire brushes or abrasive sponges and cleaning agents that may damage the product surfaces.
- Do not use detergents with solvents that can affect the structure and consistency of plastic (benzene, toluene, acetone, etc.).
- Clean the stainless steel areas of the product with maintenance spray, the pH value of the cleaning spray: 10.2.
- The density value of the maintenance spray should be 0.855 g/cm³.
- Must be biodegradable.
- The cleaning spray must not contain AOX.

15.2. Disinfection

- The disinfectant Detrosept AF applied in the disinfection process must be a fast-acting alcohol-based spraying and wiping disinfection product that does not contain phenol and aldehyde.
- Disinfecting use is applied by spraying to completely cover the pre-cleaned medical device (spray distance 30 cm). For its effectiveness, it should be kept for the duration of microbiological activity and the product should be wiped by choosing a sterile, non-particulate cloth.
- **Properties of the disinfected product;**
 - ✓ Effective in 1 minute (bactericidal, fungicidal, virucidal, tuberculocidal)
 - ✓ Aldehyde and phenol free
 - ✓ Compatible with glass, ceramic, silicone, plastic (including plexiglass), wood, aluminium and stainless steel materials
 - ✓ Broadspectrum of action
 - ✓ 10% ethyl alcohol, 20% propane 2-ol, 0,25% Didesylmethylpoly(oxyethyl) ammonium propionate, preservative additives, perfumed deionised water
 - ✓ didecylmethylpoly (oxyethyl) ammonium propionate, protective
- **Storage conditions of the disinfected product;**
 - ✓ The expiry date must be 2 years from the date of production.
 - ✓ Keep the packaging tightly closed in a well-ventilated and well-ventilated place between 0-25°.

Used disinfectant with special properties.



WARNING:

Disinfectant is flammable. Keep away from sources of ignition. Irritant Avoid contact with skin and mucous membranes. In case of contact with skin, rinse with plenty of water.

16. TROUBLESHOOTING

PROBLEM	REASON	SOLUTION
BED ANY COMMAND NO ANSWER:	1. The plug is not inserted in the socket 2. The power cable is not working. 3. The control box is not working. 4. Hand Control Not Working	1. Insert the plug into the socket. 2. Replace the power cable. 3. Send the control box for repair 4. Send the hand control to service
ENGINE DOES NOT START WHILE THE SYSTEM IS RUNNING	1. Engine malfunction 2. Hand Control Malfunction	1. The motor needs to be replaced, call the service centre. 2. Send the hand control to service
THE ENGINE DOES NOT START WHEN THE SYSTEM IS RUNNING, FROM THE CONTROL BOX " CLICK " IF THERE IS NO SOUND:	1. Control box defective 2. The remote is faulty	1. Control Box needs to be replaced, call service 2. The remote needs to be replaced. Call the service centre.

17. SAFE STORAGE OF THE PATIENT BED

To prevent damage to the patient bed during storage;

- ✓ Wrap the power cable around the cable hook on the patient bed.
 - ✓ Remove the patient bed accessories and position them horizontally on the bed.
 - ✓ Pack the patient bed so that there are no moving parts due to involuntary shaking.
 - ✓ Keep the height of the patient bed to a minimum.
 - ✓ Adjust the position of the patient bed so that it is level.
 - ✓ Switch the brakes to the off position.
 - ✓ Do not place the bed under load during storage.
 - ✓ Wrap electrical systems with protective packaging materials.
-
-

18. ENVIRONMENT

18.1 Environmental Protection

PEDIA PALS carries the idea of taking great steps in protecting the environment by thinking about future generations. The materials of this product are environmentally compatible. It does not contain hazardous substances. Attention should be paid to the symbols and signs as information is provided with the symbols and signs used. The product uses recycled steel, electrical components, packaging materials, plastics and wood materials.

The maximum acoustic noise level to the environment when the patient bed is in operation is **60 dB**.

- ✚ Defective electrical materials must be returned to the manufacturer without being discarded.
- ✚ Return the packaging material to the manufacturer or a licensed recycling company.
- ✚ Return the used defective plastic material to the licensed recycling company or the manufacturer.

NOTE: For necessary information, support should be obtained from our after-sales service unit.



19. WARRANTY

PEDIA PALS shall only be held responsible for regular servicing and product reliability.

This product is covered by a **24-month** warranty from the date of purchase. Defects caused by production and assembly defects are free of charge. Defects caused by usage errors are not covered by the warranty. Product use is determined within the framework of the terms and conditions determined by the standard.

20. CONTACT

MANUFACTURING COMPANY

BRAND: PEDIA PALS

ADDRESS: 230 Grider St. Buffalo, NY 14215 USA

PHONE NO: 1-888-733-4272

E-MAIL: sales@pediapals.com

TYPE: Hospital Furnishings And Equipment For Medical Purposes



MODEL: 600026

LIFE OF USE: 10 YEARS

SERVICE STATION: Pedia Pals