


Control Panel Guide

Product Code- FC-PHR0027



Therapy Mode - Press to toggle between Alternate, Pulsation and Static modes on the mattress.	Panel Lock - Press and hold for 3 seconds to lock/unlock all controls from patients or other non-caregivers.
Patient Care/Maximum Inflation - Sets the maximum pressure	Power Failure Alarm - Indicates loss of power to control unit.
Cycle Time - When mattress is in Alternate Mode, this feature allows you to select the duration of alternation in 3-95 minute increments. When mattress is in turning mode, this feature allows you to select the duration of turning in 3-95 minute increments.	Sensor Disconnect - Indicates the Auto Fowler sensor has been disconnected. Please check the connection and inspect sensor cable for damage.
Power Switch - After the unit is plugged in, the LED will illuminate AMBER. Press the power switch to operate the mattress, the LED will turn GREEN.	Patients BMI Setting - Press to toggle for English/Metric measurements, Patient's height setting, Patient's weight setting.
Audio Alarm ON/OFF - Press to silence alarm until bed resets or alarm issue resolved.	Comfort Setting - Manually choose from selectable levels of air pressure for support of the patient. Manufacturer recommends using a hand check to select proper pressure setting for an individual patient.
Auto Fowler - Automatically detects the head angle of the mattress and and increases pressure for patient support.	Rotation Therapy Selection - Press Turn to select left turn only, right turn only or full turn
Auto Set - Automatically sets pressure based on patients height, weight and position.	LCD Display - Shows active Control unit settings.
Manual Therapy Selection - Changes English/Metric measurements, Patients height setting, Patients weight setting.	

Operational Guidelines

- Remove the existing mattress from the bed frame.
- Replace the standard mattress with mattress replacement system (orient mattress so that the air tube is at the foot of the bed).
- Secure mattress straps attached on the bottom cover to the bed frame.
- Position the control unit on the foot board of the bed frame.
- Attach the air tube connector and Auto Fowler sensor cable to control unit's socket.
- Verify that air hoses are not kinked under the mattress.
- Check the CPR to ensure it is set to "Close" position.
- Plug in the control unit power cord located on the left side (the STANDBY LED will illuminate AMBER).
- Press the STANDBY/OPERATE switch button on the control panel (OPERATE LED will now be illuminated GREEN and the control unit will be in operation).
- CPR Deflation: The air hose connectors must be disconnected from the controller to quick release the air in an emergency when CPR is to be performed. The RED CPR valve located at Head End of bed must be twisted open to quick release the air in the Raised Side Rails when in an emergency where CPR is to be performed.

Cleaning

The mattress should be cleaned on the bed weekly using a damp soft cloth and mild detergent. If top cover or base cover becomes grossly soiled, put on clean gloves, plastic gown and eye protection before removing top and base covers and dispose of according to standard hospital procedures for contaminated waste and replace with clean covers.

Covers can be washed and thermally disinfected in a washing machine. Never use phenol based cleaning solutions.

Duet Suggested Weight Guideline

Level	*Suggested weight range (KG)	*Suggested weight range (LB)
P1	40-64kgs	88-140lbs
P2	64-87kgs	140-192lbs
P3	87-111kgs	192-244lbs
P4	111-134kgs	244-296lbs
P5	134-158kgs	296-348lbs
P6	158-181kgs	348-400lbs
P7	181-205kgs	400-452lbs
P8	205-229kgs	452-504lbs
P9	229-252kgs	504-556lbs
P10	252-272kgs	556-600lbs

**Please note that the above charts are just a guideline, not a rule. Because of different body types, shapes and height we recommend using the hand check method to ensure that the pressure setting is optimal and that the patient isn't bottoming out.*