

# REHAB WALKING TOTAL SUPPORT SLING SIZING GUIDE

(formely Prism Medical)

To find waist circumference, measure around the patient 2" above the navel.

<b>Size</b>	<b>Waist Circumference</b>
<b>Junior</b>	<b>25"-32"</b>
<b>Small</b>	<b>32"-36"</b>
<b>Medium</b>	<b>37"-52"</b>
<b>Large</b>	<b>53"-60"</b>
<b>X-Large</b>	<b>61"-68"</b>
<b>XX-Large</b>	<b>69"-76"</b>

Remember to use your best clinical judgement when choosing a sling size which best fits your patient. People come in all shapes and sizes and two people who may weigh the same or be the same height but require different sized slings.

If patient weighs over 400lbs a sling consultation is recommended.  
Please call (866) 891-6502