



*Separate chambers to assure stability*

Re-Order #	Description	Size	Unit
751046	Stability Plus Gel Foam Cushion w/Low Shear I Cover	16" x 16" x 2.5"	Each
751047	Stability Plus Gel Foam Cushion w/Low Shear I Cover	18" x 16" x 2.5"	Each
751048	Stability Plus Gel Foam Cushion w/Low Shear I Cover	20" x 16" x 2.5"	Each

### **Purpose**

This is the ultimate cushion that combines stability, pressure relief and positioning. The separate pods enhance stability while the raised pommel helps control sliding.

### **Important Safety Feature**

The Stability Plus Cushion has a cloth outer cover designed with safety straps that must be secured to the wheelchair to prevent the cushion from sliding off the seat.

### **Easy Care**

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. The Low-Shear outer cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.

*This product is guaranteed for 1 year.*