INSTRUCTIONS:

Nuprodx bath/commode chairs and slider systems are known for their ability to fit in cramped spaces others simply can't. While we have yet to find a bathroom we can't accommodate, a few simple measurements will insure the best fit of your chair. Follow the descriptions below, and enter the measurements. Round up or down to ½-inch increments, and use the decimal system (for example, 10.5 for 10 and one-half inches). If you have any questions, please call and let us help you, because that's what we do!

- A Overall length of bathtub
- B Length of bathtub floor
- **C** Width of bathtub floor (flat portion only—do not include any rounded area)
- D Overall width of bathtub
- E Distance between edge of toilet bowl and bathtub
- F Widest portion of toilet bowl (including "wedge" shape up to face of tank
- **G** Distance between edge of toilet bowl and closest object/wall
- H Distance from front edge of toilet bowl to closest object/wall

J *When the toilet is NOT next to the tub: Distance from the center of the bathtub to the edge of the outside bathtub ledge

- K Distance from center of toilet bowl to center of bathtub drain/spout
- L Distance from top edge of bathtub to floor of bathtub
- M Distance from top edge of bathtub to bathroom floor
- **N** Distance from top of toilet bowl to bathroom floor (with seat and lid in raised position)



